


*Presenting*

# JUPITOR

A Complete Range of Centrifuges  
for the Dairy Industry



MINICENTRO



HAND OPERATED



MIDICENTRO

**Leader of The New Generation Centrifuges**



... Insist On

**JUPITOR**

Because You Deserve The Best

**GERBER**  
CENTRIFUGES



Model- JSEL 11

**MINICENTRO**

- Electrically operated
- Dual Purpose : 8 Butyrometers or 4 Sol. Tubes
- Timer (0 - 5 Minutes)
- Mechanical Brake
- Protective cover
- Portable with Anti-vibration rubber pads



Model-JSHO01

- Basic Hand operated centrifuge
- Capacity : 8 / 12 / 24 Butyrometers
- Smooth running with perfect balance
- Acid proof clamps and sockets

**MIDICENTRO-EX**

- Electrically operated
- Capacity : 8 / 12 / 24 Test
- Timer (0 - 5 Minutes)
- Mechanical Brake
- Protective Cover
- Stainless Steel Disc
- Export quality



Model- JSEL 24

**MIDICENTRO**



Model-JSEL22

- Electrically operated
- Capacity : 8 / 12 / 24 Butyrometers
- Timer (0 - 5 Minutes)
- Mechanical Brake
- Protective Cover
- Portable with Anti-vibration rubber pads



Model-JSEL23

**MIDICENTRO-S**

- Electrically operated
- Capacity : 8 / 12 / 24 Butyrometers
- Timer (0 - 5 Minutes)
- Mechanical Brake
- Protective Cover
- Speed Indicator
- Special Geared Motor at 1100 rpm (ISO Standard)

