


Deepvision
PLUS
MICROSCOPES



Deepvision
PLUS
DX-10

- Sharp, bright and clear images through high quality optics
- Perfect for teachers and students alike
- Great for viewing biological specimens on slides
- High power quality compound monocular microscope



Deepvision
PLUS
DX-10



DX-10

▼ Technical Specification

Head	Monocular Head 45° Inclined 360° Rotatable
Eyepiece	WF 10x
Stage	Plane Stage With Paired Clips
Objectives	Achromatic Objectives 4x, 10x & 40x (SL) and 100x Oil (Optional)
Nosepiece	Triple Revolving ball bearing Nose Piece for life long smooth rotation
Focusing	Coarse Focusing Adjustment Mechanism on ball bearing
Condenser	Six stop Dial Diaphragm with Single Lens Condenser N.A. 0.65
Illumination	Built-in-illumination 1W LED Lamp with Rechargeable Device
Safety Device	Adjustable down stopper for preventing accidental damage to slides
Packing	full moulded thermocol container

