


Deepvision
PLUS
MICROSCOPES



Deepvision
PLUS
DM-R Series
Research Microscope



Deepvision
PLUS

DM-R Series Research Microscope



DM-RL



DM-RS

▼ Technical Specification

Body	Alluminium die cast parts with all critical movements, based on ball-bearing guides	
Observation Head	Binocular Head 45° inclined, 360° rotatable	Sidentopf Binocular Head 30° inclined, 360° rotatable
Stage	Co-axial double plate mechanical stage with ball bearing slides for frictionless and jerk free movement	
Condenser	Moveable condenser (NA 1.2) with IRIS diaphragm and filter holder on rack and pinion	
Focusing	Coaxial Coarse and Fine focusing mechanism	
Objectives	SEMIPLAN objectives 4x, 10x, 40x (SL) and 100x (SL) OIL (Anti fungus)	
Eyepiece	WF 10X (F. no. 18) (Pair)	
Nosepiece	Quadruple reverse revolving ball bearing nosepiece with rubber grip for long-life smooth rotation	
Safety Device	Adjustable down stopper for preventing accidental damage to slides	
Illumination	6V 20W Halogen bulb with variable illumination control	
Packaging	Full moulded thermocol container	

OPTIONAL ACCESSORIES

1. Trinocular Head / Sidentopf trinocular head
2. DIN Plan Achromatic objectives 4x, 10x, 40x (SL) and 100x (SL) OIL
3. LED Illumination System
4. Phase contrast attachment

Deepvision
PLUS
MICROSCOPES

