

*Deepvision*  
PLUS  
MICROSCOPES



*Deepvision*  
PLUS  
**DMM-11 Series**  
Pathological Microscope

- High power quality compound monocular microscope
- Great for viewing biological specimens on slides
- Perfect for teachers and students alike
- Sharp, bright and clear images through high quality optics



# Deepvision PLUS DMM-11 Series Pathological Microscope



**DMM-11 A**



**DMM-11 B**



**DMM-11 C**



**DMM-11 D**

▼ Technical Specification

Body	Stable and robust, assembled from alluminium die cast parts with heavy base			
Observation Head	Monocular tube inclined upto 90°	Monocular head 45° Inclined, 360° Rotatable	Monocular tube inclined upto 90°	Monocular head 45° Inclined, 360° Rotatable
Stage	Detachable mechanical stage		Co-axial mechanical stage with ball bearing slides for frictionless and jerk free movement	
Condenser	Moveable condenser (na 1.2)with IRIS diaphragm and filter holder. on rack and pinion			
Focusing	Separate Coarse & Fine focusing adjustment			
Objectives	Achromatic objectives 10x, 40x (SL) and 100x (SL) OIL		Semi-plan objectives 10x, 40x (SL) and 100x (SL) OIL (Anti fungus)	
Eyepiece	WF 10x and H5x or H15x			
Nosepiece	Triple revolving nosepiece with precise centering and positive click stops.		Ball bearing nose piece with rubber grip for life-long smooth rotation	
Safety Device	Adjustable down stopper for preventing accidental damage to slides			
Illumination	Plano concave mirror dia meter 50mm mounted on a gimble			
Packing	Full moulded thermocol container			

**OPTIONAL ACCESSORIES**

1. Sub stage pin type lamp working on 220 volts - 15 watt.
2. WOODEN BOX with lock and key.
3. Objective 4x

**Deepvision  
PLUS  
MICROSCOPES**

