

DIGITAL TABLET HARDNESS TESTER

Tablet breaking force (commonly known as Tablet Hardness) is a measure of mechanical integrity of a tablet. It is widely used in pharmaceutical industries and research institutions.

This tester is designed for determination of breaking force of tablets of different shapes quickly and accurately. Tablet breaking force is displayed on a LCD display.

SALIENT FEATURES

- * Digital Display of Breaking Force (Hardness)
- * Easy to Use Table Model
- * Suitable for up-to 30mm tablets
- * Easy to Calibrate



Description	SHT-293
Type	Table Model with Digital Display
Tablet Size	Up to 30mm length or diameter
Force Range	0 to 30 kg force (0-294 Newton)
Display	Digital (LCD)
Minimum Force Readability	0.1kg
Application of Force	With screw knob
Zero Setting	with a press of key

*Due to continuous development, look and specifications may change