



Powerful and easy-to-use RM 02 Plus mini centrifuges are small enough so that each workstation can be equipped with a personal centrifuge.

- Brushless DC (BLDC) Motor
- Unique Closed rotor reduces friction to ensure better airflow, low heat generation & whisper quiet operation
- Max speed of 6000 RPM - 2000 x g - can process upto 8 microtubes of 1.5/2 ml at a time
- Electronic brake for immediate stop on lid opening for user safety
- Save time by doing quick runs by simple opening and closing the lid
- Microprocessor controlled for greater efficiency - Next gen design with digital display
- Robust construction for long life under repeated usage
- User calibration facility
- Carbon isolation feature (dust does not enter rotor chamber)
- Replaceable strip rotor & adaptors for 0.2 ml & 0.4/0.5 ml microtubes
- Extremely compact footprint saves valuable bench space- suitable for use in fume hoods or cold rooms
- Settable speed (upto 6000 RPM) and time (upto 99 minutes)

Technical Data

Specifications	RM 02 Plus
Rotor Capacity	8 place closed rotor for 1.5/2.0ml microtubes
Speed RPM / RCF (g)	6000/2000 g
Digital Display	Yes
Electronic Safety Brake	Yes
Noise Level	< 55 dB
Power Supply	Universal (95 to 230V, 50/60Hz)
Dimensions (WxDxH) in mm	162 x 157 X 115 mm
Rubber suction pads	Yes
Rotors and Adaptors (included in standard pack)	8 Slot x 1.5/2.0ml Closed Rotor
	PCR Strip Rotor (2 x 8 x 0.2 ml)
	Reduction Adaptors for 0.2 / 0.4 ml microtubes

