

Digital Friability Test Apparatus

(BP/USP Standard)

DURALAB

Salient Features

- * Programmable Rotation Time with Digital Display
- * Auto Removal of sample without opening the drum
- * Auto stop & Alarm at completion of Test
- * Special Low rpm motor for long life



| Description | SFA-263 | SFA-264 |
|------------------------|--------------------|--------------------|
| Rotation Speed | 25+-1 rpm | 25+-1 rpm |
| Time Range | 1 Sec. To 99.59min | 1 Sec. To 99.59min |
| No. of Drum | One | Two |
| Electric Motor | Low rpm A.C. Motor | Low rpm A.C. Motor |
| Audio Alarm | Yes | Yes |
| Operation Time | Programmable | Programmable |
| Auto Removal of Sample | Yes | Yes |