



REMI Water Bath RWB Series



RWB 6

REMI Water Bath is multipurpose water bath with adjustable double line six holes to hold different size samples for easy & hygienic heating. RWB provides precise and constant temperature helps achieving uniform heating for the samples.

Salient Features:

- Microprocessor control with timing function
- Digital display
- Audible and visual alarm for over temperature

Technical Specifications

| Parameter | RWB 6 |
|--------------------------------|--|
| Temperature range | RT \pm 5 $^{\circ}$ C ~99 $^{\circ}$ C |
| Temperature stability | \pm 0.5 $^{\circ}$ C |
| Tracking alarm | \pm 2 $^{\circ}$ C |
| Chamber volume | 15 Ltrs |
| Interior Dimension (WxDxH, mm) | 450x300x110 |
| Timing range | 1~999min |
| Remark | Double lines six holes |
| Electrical requirements | 220V, 50Hz, Single phase |
| Power Consumption | 1500W |