

*Deepvision*  
PLUS  
MICROSCOPES



*Deepvision*  
PLUS  
DXL-20 Series  
Mini Co-axial Microscope



# Deepvision PLUS DXL-20 Series Mini Co-axial Microscope



**DXL-20M**



**DXL-20B**

▼ Technical Specification

Body	Aluminum die cast parts	
Observation Head	Monocular Head 45° Inclined, 360° rotatable	Binocular Head 45° Inclined, 360° rotatable
Stage	Co-axial double plate mechanical stage with ball bearing slides for frictionless and jerk free movement	
Condenser	Moveable Condenser (NA 1.2) with IRIS diaphragm and filter holder on rack and pinion	
Focussing	Coaxial Coarse and Fine focusing Mechanism	
Objectives	Achromatic objectives 4x, 10x,40x(SL)and 100x (SL)OIL	
Eyeiece	WF 10x and H5x or H15x	WF 10x and H5x pairs
Nosepiece	Quadruple ball bearing nosepiece with rubber grip for life-long smooth rotation	
Safety Device	Adjustable down stopper for preventing accidental damage to slides	
Illumination	6v-20w halogen bulb with variable illumination control	
Packaging	Full moulded thermocol container	

**OPTIONAL ACCESSORIES**

1. LED illumination with battery backup for uninterrupted research and learning.
2. WOODEN BOX with lock and key.

**Deepvision  
PLUS  
MICROSCOPES**

