

Disintergration Apparatus

Disintergration Digital Tablet Apparatus Disintergration Apparatus

901 & 904



1 Year
Warranty

- Digital timer with 1 second resolution
- Accurate time out provided by microprocessor based time totalizer
- Basket move 'up & Down' as per IP/USP test specifications
- Temperature controlled by sensitive thermostats

Digital Tablet Disintegration Apparatus is designed for accurate estimation of disintegration time of tablets as per IP/USP standards. The instrument 901 is designed to test two batches of 6 tablets, simultaneously & Model 904 tests single batch of 6 tablets. The unit consists of beakers having capacity of 1000 cc each. In these beakers, suitable disintegration media can be taken. These are temperature controlled using an in-built heater and a highly sensitive thermostat fitted inside the instrument. Microprocessor based time totalizer is provided for accurate time read out.

